Korean

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English

Disposable Medical Devices
Prohibition of Reuse

A Fixture Manual

■ Product information

An IBS Implant System (Fixture) functions as artificial dental root for the missing teeth of a patient suffering from a loss of a single tooth or multiple teeth. And IBS provides the patient with an opportunity to restore his/ her missing teeth by mounting abutments on them following the osseointegration with their alveolar bones. Fixture is made of medical titanium alloy (Ti-6Al-4V ELI). Regarding detailed information, please refer to manual, catalog, and/or homepage (www.ibsimplant.com). Please also refer to attached label for products code, spec, date of manufacture and/or validation.

Dental implant fixtures with a diameter of less than 3.5 mm should be used only for the replacement of the lower central and lateral incisors and the maxillary lateral incisor.

■ Indications

The IBS Implant system is intended for use in dental implant surgery. It will be surgically placed in the maxillary or mandibular alveolar bone to replace a missing natural root of tooth.

■ Sterilization

Fixture, closing screw and healing abutment have been cleaned and gamma sterilized. All other products or tools, which are not sterilized, should be steam sterilized at 132℃ for 8 minutes. It is also advised to wait 30 minutes at least after having steam sterilization. Do not use a product if packing is damaged or expired. Re-sterilization is not permitted in any circumstances and do dispose of it properly.

■ Storage conditions

Place the product in a cool, dry place, away from direct sunlight.

■ General precautions

As operational techniques for dental implant surgery require technical expertise, complex and elaborate procedures clinicians have to receive implant-related formal educations and training.

■ Preparation before use

Operating surgeon practitioner should make decisions to do implant surgery by considering bone conditions or proper surgery methods. Operating surgeon practitioner must fully check up any stomatologic, anatomical structure which possibly influences on surgery results. And they should prepare proper implant considering matters and/or factors that require attention. CT, Panoroma and/or various examinations require to exam for basic information such as antagonist and bone density and so on. Also they should establish surgery plans and make a decision appropriate fixture.

■ Precautions for use

Attempts should be made to minimize damage to the patients cellular tissue and special attentions have to be paid to temperature, potential surgical disorders, and removal of the causes of contaminations and infection. Implant surgery requires higher degrees of accuracy and attentions. A deviation form the principles of implant surgery will increase a danger of failure in osseointegration. A space for placement of a fixture agreeable to the patient's bone quality and fit to the fixtures specications must be secured. Special care should be taken for fear that the connection part may have trouble in disconnection or get damaged it a fixture is placed beyond 50Nem torque. Be sure to tighten the closing screw with a hand torque 5~8 Nem) after the fixture is completely settled: connection with a higher than specified torque may not permit smooth separation. An implant placed at an angled of 30° or higher may result in fracture. Pay attention to small-size implant component and tools as they may be swallowed by the patient during surgery. The bone quality after the placement of a fixture and its initial fixation play an important role in the decision on loading time.

Radiograph can be used to observe and diagnose the patients conditions including bone loss around the implant and veridical changes in osseointegration. It the observation finds that the imolant moves or there is a bone loss of more than 50%, a judgment must be made on whether to remove the implant.

■ Precautions for patients

It is important to practice good oral hygiene. Avoid imposing too much pressure until the prosthesis is finally mounted. Be sure to follow and abide by the instructions and guidelines given by yours dentist.

■ Warnings

IBS implants are medical devices. Do not use them for any other purpose than originally intended: do not let anyone who is not authorized handle them: do not modify or alter them in any manner. Any person who uses them must receive special training sessions so as to use them in safe and effective manner for technologies necessary for dental implant procedures are surgically specialized and requires complex process. Inaccurate diagnosis on a patient and wrong selection of a treatment lead to tailure in implant surgery and result in bone loss. Use of an electric surgical appliance or a surgical laser around a metallic dental implant is prohibited for fear of electric shocks or burns. Loose implants, bone loss and chronic infection can end up in failure

in implant surgerv.

No allowed angular correction of the abutments nor placing the implant bodies at an angle (i.e., must be used in parallel to the direction of occlusal loading forces)

■ Side effects

Implant surgery may later be subject to several problems including loss in primary stability of the implant and loss of the prosthesis. The loss of primary stability can result from such factors as quantitative and qualitative deficiency in remaining bone, patient's poor oral hygiene or uncooperative altitude, implant's movement, deterioration of local tissue, and incorrect position and alignment of implant.

■ Contraindications

Patients who have any of the following conditions must not receive a surgical treatment:
Oral contradictions: pathological alveolar ridge, dry mouth, previous radiation therapy on the jawbone, pachy glossia, pathological changes of oral mucosa (Vitiligo, lichen planus, or stomatitis) and poor oral hygiene
Contraindications for temporary prohibitions of the surgery: when a patient takes specially prescribed medicines such as anticoagulant and immunosuppressant, has been under physical and mental stress, has contracted acute inflammatory or airborne diseases, or is pregnant.

Psychological and social contraindications: when a patient is uncooperative, is addicted to alcohols and/or drugs, suffers from neurosis, or has been hospitalized in a mental clinic.
Medical contradictions: when a patient takes prescribed medications use of steroid drugs for long-term use of antibiotics, etc) or suffers from diabetes, cardiac disease, metabolic disorder such as congestive heart failure, and blood disorder. And patient who is judged not to be in a conditions for surgical procedure.

*MR Statement

The system have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the System in the MR environment is unknown. Scanning a patient who has one of these devices may result in patient injury.

■ How to use

Follow general guidelines for implant surgery and refer to our catalogue and surgery manual for more information.

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Innobiosurg, LTD. 44-19 Techno 10-ro, Yuseong-gu, Daejeon, Korea